


## Offset Acetabular Cup Impactor



Our offset impactor provides controlled positioning and impaction of acetabular cups through a small incision, enabling a less invasive surgical approach.

### DESIGNED TO PERFORM

- Impaction load applied in-line to acetabular cup.
- Remote cup orientation.
- Unique locking mechanism for accurate cup placement during impaction; the cup will not rotate once locked in place.
- Integrates Greatbatch Medical's signature EZ Clean features.

### CUSTOMIZED TO YOUR NEEDS

- Customizable to fit the thread profile, diameter or chamfer of your cup.
- Variants available for cups without an apical hole.
- Available with alignment guide.

Patent Pending

**Greatbatch Medical Headquarters**  
10000 Wehrle Drive  
Clarence, NY 14031  
U.S.A.  
tel +1 716-759-5600  
fax +1 716-759-2562

**US Customer Service**  
4532 Park 30 Drive  
Columbia City, IN 46725  
U.S.A.  
tel +1 260-244-6300  
fax +1 260-244-6345

**Europe & Asia Customer Service**  
L'Echelette 7  
2534 Orvin  
Switzerland  
tel +41 32 358 01 11  
fax +41 32 358 01 00

orthopaedicsna@greatbatchmedical.com

orthopaedicseurope@greatbatchmedical.com

[www.greatbatchmedical.com](http://www.greatbatchmedical.com)