




## **Q<sub>CAPS</sub>™ High Voltage Capacitors**



Q<sub>CAPS</sub> high voltage capacitors provide fast, consistent charge times for tachycardia therapy, the first time and every time. When paired with Q<sub>HR</sub> batteries, Q<sub>CAPS</sub> capacitors provide the smallest, longest-lived, highest-energy power solutions for today's tachycardia devices.

### **FAST, CONSISTENT CHARGE TIMES**

- First time, every time
- configurations 600-1000V & 30-45J
- Reduces battery load variability

### **OUTSTANDING EFFICIENCY**

- Highly energy efficient over life
- Up to 30% smaller than aluminum electrolytic
- No reform required

Q power solutions from Greatbatch Medical combine fully characterized technology, robust design principles, world-class manufacturing and 40 years of experience to ensure unparalleled performance and reliability.

---

**Greatbatch Medical**  
10000 Wehrle Drive  
Clarence, NY 14031  
tel 716-759-5600  
fax 716-759-2562

**sales@greatbatchmedical.com**  
**www.greatbatchmedical.com**